



Boise Pathology Group, PA ~ 190 East Bannock Street ~ Boise, Idaho 83712
Phone: (208) 381-2367 ~ Fax: (208) 381-4762
www.boisepathologygroup.com

SPECIMENS **NEVER** TO BE PLACED IN ROUTINE FORMALIN FIXATIVE

Certain specialized procedures require that the specimen not be fixed in formalin. These include:

- A. Cytogenetic studies
- B. Gout Crystals
- C. Flow Cytometry
- D. MOH's
(Contact Caroline Durupt @ (208) 761-0363 to schedule)
- E. Frozen section analysis
(Contact St. Luke's histology @ (208) 381-2370 to schedule.)
- F. Immunofluorescence studies
- G. Nerve/Muscle biopsies
- H. Cytology fluids

Please call St Luke's histology at 208-381-2370 with questions.